

## Features

- ▶ Integrates with CE Series drop-in ceiling panel.
- ▶ Accessible control and service isolation.
- ▶ Multi panel control from one location.
- ▶ Process filtration, standalone POU and service lockout options for core or clean (cGMP) side.
- ▶ Customizable panel configurations.
- ▶ Reduce on-site cost, labor, and risk.

## Technical Data

### Temperature Rating:

50°F (10°C) to 140°F (60°C)

### Pressure Rating:

500 psig (35 bar) inlet

### Regulator Outlet Rating:

0-150 psig (0-10 bar)

### Inlet/Outlet Size:

1/2"

## Ordering Information

### UPSS-4A40

#### Gas Connection (Qty/Service)


- A - Compressed Air
- N - N<sub>2</sub>
- O - O<sub>2</sub>
- C - CO<sub>2</sub>

#### Panel Material


- UPSS - 304 Stainless Steel (brushed)
- UPCS - Carbon Steel (powder coated)

## UP Series Enclosure

### Multi-use Inert Gas Wall Mount Enclosure



## Dimensions



Front

Side

**Note** Product availability and specifications contained herein are subject to change without notice. Consult factory for potential revisions.



**BioPharm Engineered Systems**

A Critical Process Systems Group Company

200 Bulfinch Drive, Andover, MA 01810

+1 978.691.2737

F +1 978.691.2733

| [www.bpesys.com](http://www.bpesys.com)

**CONFIDENTIAL**

REV: CPS 03032016