

Features

- ▶ Integrates with CE Series drop-in ceiling panel.
- ▶ Accessible control and service isolation.
- ▶ Multi panel control from one location.
- ▶ Process filtration, standalone POU and service lockout options for Core or Clean (cGMP) side.
- ▶ Customizable panel configurations.
- ▶ Reduce on-site cost, labor, and risk.

Technical Data

Temperature Rating:

50 to 140°F (10 to 60°C)

Pressure Rating:

500 psig (35 bar) Inlet

Regulator Outlet Rating:

0-150 psig (0-10 bar)

Inlet/Outlet Size:

1/2"

Ordering Information

UPSS-4A4N4O4C


Gas Connection (Qty/Service)

- A - Compressed Air
- N - N₂
- O - O₂
- C - CO₂




Panel Material

- UPSS - 304 Stainless Steel (Brushed)
- UPAL - Aluminum (Powder Coated)

UP Series Panel Multi-use Inert Gas Wall Panel



Dimensions



Note Product availability and specifications contained herein are subject to change without notice. Consult factory for potential revisions.



BioPharm Engineered Systems

A Critical Process Systems Group Company

200 Bulfinch Drive, Andover, MA 01810

+1 978.691.2737

F +1 978.691.2733

| www.bpesys.com

CONFIDENTIAL

REV: CPS 03032016