

F1 Series

Advanced Technology Division

Features

- ▶ Fully customizable to meet requirements.
- ▶ Modularized design to streamline on-site construction.
- ▶ Critical service cleanroom fabrication.
- ▶ NEHP Engineering, Design, and Logistical Services are available to support through the construction and commissioning phases.
- ▶ Reduce on-site cost, labor, and risk.

Technical Data

Frame Material:

304L Stainless Steel / Aluminum

Frame Finish:


#4 Brushed or Powder coated White (std),
Color match available

Pipe/Tube Specs:


Process dependent

Pipe/Tube Materials:

PVC, PP, PVDF, PFA, CS, SST standard,
other materials available



Dimensions



Ordering Information

Consult factory with details to discuss ordering information.

Note Product availability and specifications contained herein are subject to change without notice. Consult factory for potential revisions.

 **BioPharm Engineered Systems**
 A Critical Process Systems Group Company
 200 Bulfinch Drive, Andover, MA 01810

+1 978.691.2737 F +1 978.691.2733 | www.bpesys.com

CONFIDENTIAL

REV: CPS 03032016